

TO WHOM IT MAY CONCERN

**Date** 29/01/2019  
**Hospital name** NH Hospital a.s. - Nemocnice Hořovice  
**Street** K Nemocnici 1106/14  
**City** Hořovice, ZIP:268 01  
**Country** Czech Republic  
**Phone** +420 775 324 006  
**Contact person** Ing. Kamil Kudelka  
**Title** Mr.  
**FSN monitor model** FS-S4601DT, FS-P2607D

We hereby confirm that in period of 2017-2018, the Hospital resp. Central operation rooms acquired FSN Monitor models: 1 x FS-S4601DT, 2x FS-P2607D, manufactured by Foreseeson.

These monitors are connected via DVI-D to Olympus OTV-S190 and ViseraELITE II HD2D/3D videoprocessors which and is used in the daily clinical practice of our institution when working with laparoscopes and endoscopes.

On average, every day during the various procedures using Monitors for average **7-9 hours a day**, average of **10 patients** are examined **daily**.

**More comments**

FS-P2607D are used mostly for gynecology and urology operations (3D) as looped two monitors, one installed on right side on the tower and one the left side on the ceiling arm. Connected together by DVI-D cable (ClickTronic HQ OFC, DVI-D 24+1 - DVI-D 24+1, M/M, Dual Link, Au connectors, 20m)

We did not have any technical issue so far. Brightness, colors are stable as I would expect. 3D works well.

FS-P2607D was tested against SonyLMD-X310ST before the purchase and we did not see such a big different except price and screen size, so we decided for FS-P2607D and according to doctors they are satisfied enough.

FS-S4601DT is used not only as a third monitor for student watching long and complicated operations like IBD operation is but it serves as well as gate to allow the PACS for x-rays, CT, MRI, US images consultation and review of previous operations right in the room.

It is connected to the 3G SDI out of FS-P2603D and by cable Creation Bi-Directional HDMI Female to DVI-D(24+1) Male Adapter, 1080P DVI to HDMI Converter, 3D to the PC stick (PACS access).

It works well. We did not have yet any technical issue as well.

**During the usage of the monitors there were no complaints and negative reviews from medical workers.**

**Contact person** Ing. Kamil Kudelka  
**Title** IT Director  
**Signature**

**Stamp**  
NH Hospital a.s.  
Okružová 1135/ 44  
155 00 Praha 5  
IČ: 278 72 963 ①



**Date** 29/1/2019